## Alaska Scientific Crime Detection Laboratory

## DataMaster Repair Form Effective: 10/10/2022 Version: 1.0 Analyst(s): Carlie Glaister Instrument #: 10068 Start Date: 1/130/23 Reason for Repair Internal standard errors Troubleshooting/Repair Cleaned the internal standard with 184. ✓ Verify software version VOC (attached) 100 Internal Standard value = 3.00% Non-Drinking Subject Test (attached) 15-14.16 1. Not performed due to failed VOC. Ready for Adjustment Ready for Calibration Further Repair Needed Return to Factory Proceed to Close Out End Date: 11/30/23

Page 1 of 1

All printed copies are uncontrolled.

Approved by: Chemistry Supervisor

## VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Serial #: 100687 Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Date: 11/30/2023

External Standard Test Values	Diagnostic Check	
	VERSIONS	
	DMT: 3.02	
	PIC: 3.02	
	Modem: 2.6	
	Questions: 2.2	
	TEMPERATURES	
	Sample Chamber = 48.8°C	PASSED
	(44.0°C - 52.0°C)	
	Breath Tube = 46.5°C (38.0°C - 50.0°C)	PASSED
	SETTINGS	
	Lamp Voltage = 1.56 V (0.94 V - 2.18 V)	PASSED
	Cooler Voltage = 1.56 V	PASSED
	(0.94 V - 2.18 V) Bias Voltage = 80 V	PASSED
	(48 V - 112 V)	PASSED
	Chopper Freq = 534 Hz	PASSED
	(475 Hz - 575 Hz)	
	Barometer = 28.8 in	
	PUMP INFO	
	Flow Rate = 4.403 L/M	PASSED
	(3.500 L/M - 6.500 L/M)	
	DETECTOR INFO	
	PUMP ON	PASSED
	(0.002 V <= 0.010 V)	
	PUMP OFF	PASSED
	(0.002 V <= 0.010 V)	
	FILTER INFO	
	Filter 1	PASSED
	Filter 2	PASSED
	Filter 3	PASSED
	INTERNAL STANDARD	FAILED-
	Xq = 0.097 14.16% (0.00% - 0	

Diagnostic Check Failed