Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form

Version: DMEF2015 R0

E.	fective: 7/10/2015		Status: Active
Instrui	ment # 100416		
Date:	9/14/15	Analyst: Brandi	Barnett
 2. 3. 	Reason for instrument evaluation Adjustment error when 1 Xq 0.0069 0.0500 <= Xq 0.0343 Visual Inspection Troubleshooting/Repair	hying to calibra	re
	Can't perform voc, detecto unto all fail.	r, filter + interna	(1 Standard (CB 9/14/15
4	VOC (attached) (1/A		
	VOC (attached) WA Internal Standard value ≤ 3.00% Yes □	No I NIA	
	Verify software version		
	Return to Factory for Repair		
	Ready For Calibration		
	Ready For Certification		

Issued: 7/10/2015

CALIBRATION FACTORS

Alaska Department of Public Safety

DATAMASTER dmt: 100416

Date: 08/12/2015 Time: 10:15:00

OPERATOR NAME: COLLEEN S OBRYANT CERTIFICATE #: 6237

Current Cal factors

CALIBRATION FACTORS

Alaska Department of Public Safety

DATAMASTER dmt: 100416

Date: 09/14/2015 Time: 14:23:03

OPERATOR NAME: BRANDI BARNETT

CERTIFICATE #: 7038

Ca = 0.1000

Adjustment Error

Aftempted calibration

CALIBRATION FACTORS

Alaska Department of Public Safety

DATAMASTER dmt: 100416

Date: 09/14/2015 Time: 14:25:26

OPERATOR NAME: BRANDI BARNETT CERTIFICATE #: 7038

Adjustment Error

Attempted calibration

VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Date:09/17/2015

External Standard Test Values	Diagnostic Check		
	VERSIONS DMT: 3.01 PIC: 3.02 Modem: 2.6 Questions: 2.2 TEMPERATURES Sample Chamber = 48.9°C PASSED Breath Tube = 44.0°C PASSED PUMP INFO Flow Rate = 4.557 L/M PASSED DETECTOR INFO PUMP ON PUMP OFF FILTER INFO Filter 1 Filter 2 Filter 3 INTERNAL STANDARD Diagnostic Check Failed		

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

Serial #: 100416

(1) I am a Forensie Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Nita J.Bolz Scientific Director State Breath Alcohol Program			
Subscribed and sworn before me this	day of	, 20	

Technician Screen

DATAMASTER dmt: 100416

Date: 09/16/2015 Time: 3:21 PM

Temperatures (°C)

Sample Cell: 48.8°C Breath Tube: 46.7°C

Current Barometer 29.2 in Volume (Ltr) 0.00

Settings

Lamp: 1.60 V
Bias: 82 V
Cooler: 1.64 V
Chopper: 533 Hz

Voltages (V)

Flow:

Detector:

1.07-V 0.356 V

Voltage @ Alter 1 + quarta

Technician Screen

DATAMASTER dmt: 100416

Date: 09/16/2015 Time: 3:21 PM

Temperatures (°C)

Sample Cell: 48.8°C Breath Tube: 46.7°C

Current Barometer 29.2 in Volume (Ltr) 0.00

Settings

Lamp: 1.60 V
Bias: 82 V
Cooler: 1.64 V
Chopper: 533 Hz

Voltages (V)

Flow: Detector: 0.154 V

voltage Ofilter 1