Alaska Scientific Crime Detection Laboratory

DataMas	ter Evaluation Form
Issued: 05/15/2020 Effective: 05/15/2020	Version: DMEF 2020 R0 Status: Active
Instrument #	
Start Date: <u>07/13/21</u>	Analysts: <u>Deren Watton</u>
Reason for instrument evaluation	
Filter wheel errors	
Visual Inspection Clean	
	Replace missing screws
Cermetex Box is secure	Ensure modem cable present
Gas Valve upright or secured to wal	All tubing and connections are secure
O-ring inspection and ensure 2 O-ring	ngs in replacement bag
Troubleshooting/Repair	
	refurbished. Replaced oring.
Fifter wheel assembly for	refurbished. Replaced oring.
assembly maintenance,	cotory for refurbishment and filter wheel
Ran fitter test - 35 (cycles with 1 error.
VOC (attached) ☐ Internal Standard value ≤ 3.00% NA Verify software version	- filter wheel error
Return to Factory for Repair	Ready for Calibration
Ready for Certification	Proceed to Close Out
End Date: 07/13/21	
	Page 1 of 1
All printed copies are uncontrolled.	Approved by Chemistry Supervisor

VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Date:07/13/2021

External Standard Test Values

Serial #: 100363

Diagnostic Check

VERSIONS DMT: 3.02 PIC: 3.02 Modem: 2.6 Questions: 2.2

TEMPERATURES

Sample Chamber = 48.7 °C PASSED Breath Tube = 48.1 °C PASSED

PUMP INFO

Flow Rate = 4.559 L/M PASSED

DETECTOR INFO

PUMP ON PASSED PUMP OFF PASSED

Filter Wheel Error