## Alaska Scientific Crime Detection Laboratory

Data	Master Evaluation Form
Issued: 05/15/2020 Effective: 05/15/2020	Version: DMEF 2020 R0 Status: Active
Instrument #	
Start Date: 05/07/702(	Analysts: Derch walten / Carlie Bailey
Reason for instrument evaluation	
Train	1/9
Visual Inspection	
Clean	Replace missing screws
Cermetex Box is secure	Ensure modem cable present
Gas Valve upright or secured to	
O-ring inspection and ensure 2	O-rings in replacement bag
Troubleshooting/Repair	
- evaluation perform	embled for training purposes. Reassembled and ned afterward.
internal Standard	cors. No futher trunble shooting at this time
VOC (attached)	
☐ Internal Standard value ≤ 3.00% ∧	/0
Verify software version	
Return to Factory for Repair	Ready for Calibration
Ready for Certification	Proceed to Close Out
End Date: $0 \le   2  2 $	
	Page 1 of 1
All printed copies are uncontrolled.	Approved by Chemistry Supervisor

## VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Serial #: 100416 Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

External Standard Test Values

Date: 05/12/2021

Diagnostic Check		
VERSIONS DMT: 3.02 PIC: 3.02 Modem: 2.6 Questions: 2.2		
TEMPERATURES		
Sample Chamber = 48.7°C (44.0°C - 52.0°C)	PASSED	
Breath Tube = 46.5°C (38.0°C - 50.0°C)	PASSED	
SETTINGS Lamp Voltage = 1.58 V (0.96 V - 2.24 V)	PASSED	
Cooler Voltage = 1.64 V (0.98 V - 2.30 V)	PASSED	
Bias Voltage = 82 V (49 V - 115 V)	PASSED	
Chopper Freq = 531 Hz (475 Hz - 575 Hz) Barometer = 29.9 in	PASSED	
PUMP INFO Flow Rate = 4.950 L/M (3.500 L/M - 6.500 L/M)	PASSED	
DETECTOR INFO PUMP ON (0.002 V <= 0.010 V)	PASSED	
PUMP OFF (0.002 V <= 0.010 V)	PASSED	
FILTER INFO Filter 1 Filter 2 Filter 3	PASSED PASSED PASSED	
	FAILED	

Diagnostic Check Failed