VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Serial #: 100348 Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program Date:03/31/2012

External Standard Test Values

EXTERNAL STANDARD INFORMATION

NOMINAL: 0.080

Average = 0.0782Std Dev = 0.0004

TARGET AT 29.53: 0.078 LOT #: 02511080A3 EXPIRATION: 02/01/2013

BLANK TEST	0.000	18:06
INTERNAL STANDARD	VERIFIED	18:06
EXTERNAL STANDARD	0.079 1	18:06
BLANK TEST	0.000	18:07
EXTERNAL STANDARD	0.078	18:07
BLANK TEST	0.000	18:08
EXTERNAL STANDARD	0.078	18:08
BLANK TEST	0.000	18:09
EXTERNAL STANDARD	0.078	18:09
BLANK TEST	0.000	18:10
EXTERNAL STANDARD	0.078	18:10
BLANK TEST	0.000	18:11

Diagnostic Check

VERSIONS DMT: 2.00 PIC: 3.02 Modem: 2.3 Ouestions: 2.2

Questions. 2.2	
TEMPERATURES	
Sample Chamber = 48.9°C Breath Tube = 42.4°C	PASSED PASSED
PUMP INFO Flow Rate = 4.726 L/M	PASSED
DETECTOR INFO	
PUMP ON	PASSED
PUMP OFF	PASSED
FILTER INFO	
Filter 1	PASSED
Filter 2	PASSED
Filter 3	PASSED
INTERNAL STANDARD	PASSED

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Nita J.Bolz

Scientific Director

State Breath Alcohol Program

Nikki Roth, Notary Public My Commission Expires With Office



