## VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument State of Alaska

MAY 1 5 2009

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

	DataMaster cdm S/N 130579
Supe	Name Deck Sitz ID# 4094 Date S.14.09
A	Agency Anchorase Police Department Phone #
	Instrument Location Anchors Instrument Location
В	Alco S/N X 12469 Target Value . 031 High Pressure 1180
	Alco Test Values $\frac{.074}{1^{st} \text{ Alco}}$ $\frac{.075}{2^{nd} \text{ Alco}}$
	Signature Deul San Cas
	(OVER)

(Do Not write in the area below)

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.

(5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the

verification.

Scientific Director

State Breath Alcohol Program

Subscribed and sworn before the this

(Notary Seal Stamp)

Notary Public, State of Alaska Commission Expires with Office



## VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument State of Alaska

	DataMaster cdm S/N 130579
(CONTINUED I	FROM FRONT PAGE)
TAPE THE TEST STRIP FROM THE CHECK IN TH	NONDRINKING SUBJECT TEST AND THE DIAGNO HE MARKED BOXES.
NONDRINKING SUBJECT TEST	DIAGNOSTIC CHECK
ALASKA DEPARTMENT OF PUBLIC SAFETY DATAMASTER cdm 130579	ALASKA DEPARTMENT OF PUBLIC SAFETY DATAMASTER cdm 130579
MAY 14, 2009 (6) 109	MAY 14, 2009 TIME 00:37
OPERATOR'S NAME:	DIAGNOSTIC CHECK
SITZ, DEREK W OPERATOR'S NUMBER: 4096 SUBJECT'S LAST NAME:	COMPUTER: OKAY 05
VOC 1465 SUBJECT'S FIRST NAME/MI : 5 2	PROGRAM: OKAY 1415
VOC 5/14.09 O.L. #: 1234567	SOFTWARE DATE: 12/11/08
DEPT/AGENCY: ANC1 CASE/REPORT: 00-1234 TEST TYPE: V	HEATERS SAMPLE CHAMBER: 50C
ALCO TARGET VALUE: .081 ALCO S/N: X124669	BREATH TUBE: 42C BAROMETER: 30.01 in
BREATH ANALYSIS	FLOW DETECTOR: OKAY
.081 ADJUSTED FOR 30.01 in ALCO TARGET .081 00:38	PUMP
ALCO TARGET .081 00:38 BLANK TEST .000 00:39 INTERNAL STANDARD VERIFIED 00:39	HIGH SPEED: OKAY  DETECTOR: OKAY
ALCO TV 30.01 in .074 00:40 BLANK TEST .000 00:40	DETECTOR
SUBJECT SAMPLE	QUARTZ STANDARD: OKAY
ALCO TV 30.03 in .075 00:42 BLANK TEST .000 00:43	CALIBRATION: OKAY
	PRINTER TEST !"#\$\$&\()*+,/0123456789:;<=>?@ABCDEFG HIJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmno pqrstuvwxyz{ }^\(\)
1	1
	<b>↑</b>
L. C. American	