

VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm Breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

JUL 13 2009

| | |
|---|--------------------------------------|
| Supervisor/Operator Performing the Verification Procedure: | |
| DataMaster cdm S/N <u>130579</u> ✓ | |
| Name <u>Derek Site</u> | ID# <u>4096</u> Date <u>7-9-09</u> ✓ |
| A Agency <u>Anchorage Police Department</u> Phone # <u>786-8900</u> | |
| Instrument Location <u>Anchorage Jail</u> ✓ | |
| B Alco S/N <u>X124669</u> ✓ Target Value <u>.081</u> ✓ High Pressure <u>400 psi</u> | |
| Alco Test Values <u>.081</u> ✓ <u>.081</u> ✓ 1 st Alco 2 nd Alco | |
| Signature <u>[Signature]</u> ✓ | |
| (OVER) | |

(Do Not write in the area below)


I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

- (1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.
- (2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.
- (3) I am the Scientific Director of the State Breath Alcohol Program.
- (4) In that capacity, I am responsible for overseeing the breath alcohol program, including assuring that persons responsible for verifying the calibration of instruments are properly trained and qualified. I also am responsible for maintaining the records of the program.
- (5) The attached verification is a true and accurate verification of calibration that reflects a regularly conducted and regularly recorded activity of the breath alcohol program performed by a person trained and qualified to conduct the verification.

[Signature]
Nita J. Bolz
Scientific Director
State Breath Alcohol Program

Subscribed and sworn before me this 21st day of Jul 2009.

[Signature] (Notary Seal Stamp)
Carolyn M. Noland
Notary Public, State of Alaska
Commission Expires with Office



VERIFICATION OF CALIBRATION REPORT

of DataMaster cdm breath Test Instrument

State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

DataMaster cdm S/N

130579 ✓

(CONTINUED FROM FRONT PAGE)

C

TAPE THE TEST STRIP FROM THE NONDRINKING SUBJECT TEST AND THE DIAGNOSTIC CHECK IN THE MARKED BOXES.

NONDRINKING SUBJECT TEST

DIAGNOSTIC CHECK

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130579

JULY 09, 2009 ✓

OPERATOR'S NAME:

SITZ, DEREK W

OPERATOR'S NUMBER: 4096

SUBJECT'S LAST NAME:

VOC

SUBJECT'S FIRST NAME/MI :

VOC

O.L. #: 1234567

DEPT/AGENCY: ANC1

CASE/REPORT: 00-1234

TEST TYPE: V

ALCO TARGET VALUE: .081 ✓

ALCO S/N: X124669

--- BREATH ANALYSIS ---

.081 ADJUSTED FOR 30.19 in

| | | |
|-------------------|----------|-------|
| ALCO TARGET | .081 | 00:36 |
| BLANK TEST | .000 | 00:36 |
| INTERNAL STANDARD | VERIFIED | 00:36 |
| ALCO TV 30.19 in | .081 ✓ | 00:37 |
| BLANK TEST | .000 | 00:38 |
| SUBJECT SAMPLE | .000 ✓ | 00:38 |
| BLANK TEST | .000 | 00:39 |
| ALCO TV 30.19 in | .081 | 00:39 |
| BLANK TEST | .000 | 00:40 |

ALASKA DEPARTMENT OF PUBLIC SAFETY
DATAMASTER cdm 130579

JULY 09, 2009 ✓

TIME 00:34

--- DIAGNOSTIC CHECK ---

COMPUTER: OKAY

PROGRAM: OKAY

SOFTWARE DATE: 12/11/08

HEATERS

SAMPLE CHAMBER: 49c

BREATH TUBE: 42c

BAROMETER: 30.19 in

FLOW DETECTOR: OKAY

PUMP

HIGH SPEED: OKAY

DETECTOR: OKAY

FILTERS: OKAY

QUARTZ STANDARD: OKAY

CALIBRATION: OKAY

PRINTER TEST

!"#\$%&'()*+,-./0123456789;<=>?@ABCDEFGH
IJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqrstuvwxyz{|}~Δ