Alaska Scientific Crime Detection Laboratory DataMaster Evaluation Form

DataMast	er Evaluation Form
Issued: 05/15/2020 Effective: 05/15/2020	Version: DMEF 2020 R0 Status: Active
Instrument #	
Start Date: 09/13/21	Analysts: Deren Walton
Reason for instrument evaluation	
Visual Inspection	
Clean	Replace missing screws
Cermetex Box is secure	Ensure modem cable present
Gas Valve upright or secured to wall	All tubing and connections are secure
O-ring inspection and ensure 2 O-ring	gs in replacement bag
Troubleshooting/Repair 2nstrument has not been Routed filter was versum	returbished. Replaced o-ring. tely. Reloaded 3,02 software (full version)
	· · · · · · · · · · · · · · · · · · ·
Screen is very unrespons	ive and may need to be replaced. c for this reason.
Land not institute a Vo	C for This read on.
UVOC (attached) Not performed. ☐ Internal Standard value ≤ 3.00% ☐ Verify software version	DW 09/13/21
Return to Factory for Repair	☐ Ready for Calibration
Ready for Certification	Proceed to Close Out
End Date: 09/19/21	
All printed copies are uncontrolled.	Page 1 of 1
an printed copies are uncontrolled.	Approved by Chemistry Supervisor