## Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form Version: DMCIF 2017 R0 Issued: 1/13/2017 Status: Active Effective: 1/17/2017 Received From AST - St Man Instrument S/N 100418 Date 5/26/1-1. Reason for instrument return: of range high on VOC in field. 2. Change instrument and case location in JusticeTrax Breath Alcohol Database 3. Ensure change in instrument status form is in JusticeTrax 4. Upload 🔽 5. Update agency information and polling status in DM Host Calling Program 6. Check for VOCs [ **From Intoximeters** 7. Intoximeters Repair Paperwork (attached)