## Alaska Scientific Crime Detection Laboratory

	DataMaster Check	In Form
	ssued: 1/13/2017 Effective: 1/17/2017	Version: DMCIF 2017 R0 Status: Active
nstri	rument S/N_[00352	Received From Intoximeters
Date_	Ce/8/2021	Analyst Carlle Bailey
		Dewl & Watter
1.	Reason for instrument return: Returned from refurbishment. Communication issues	Instrument was howing
	communication issues.	
2.	Change instrument and case location in JusticeTrax	Breath Alcohol Database
3.	Ensure change in instrument status form is in Justice	Trax 🖸 NA
4.	Upload  NA	
5.	Update agency information and polling status in DM	Host Calling Program   NA
6.	Check for VOCs □ N 🖟	
	From Intoximeters	
7	Intoximators Panair Panaryork (attached)	

## Intoximeters

Experience · Service · Integrity

## Service Record

SRO Number: SRO-046012

Received Date: 05/17/2021

Completion Date: 05/25/2021

SRO Type: REPAIR

SRO Description: DMT Preventive Maintenance

Customer ID: C000003651

Customer PO: 050521-Walton

Contact: Derek Walton

Email: derek.walton@alaska.gov

Ship To:

Units on SRO

Alaska Department of Public Safety Procurement & Supply Section 4805 Dr. MLK Jr. Avenue Anchorage AK 99507

LEGDMT LEGACY DMT

Service Line:	**All Instruments are Calibrated to Factory Specifications**					
Unit Type:	(LEGDMT)	LEGACY DMT				
Serial Number:	50100352					
Reason:	n: Preventive Maintenance					
	planned	nned				
Reason Notes:	DMT Preventive Maintenance					
Resolution:	Prev M					
	planned					
Resolution Notes:	Removed filter wheel assembly for cleaning and lubrication. Cleaned & adjusted filter wheel solenoid. Updated filter wheel solenoid mounting screws. Updated stepper motor adjustment lock down screw. Replaced filter wheel motor pinion. Cleaned, lubricated, & adjusted filter wheel moving parts. Applied epoxy to stepper motor wires. Removed, inspected, and applied metallic seal tape to n blocks of Sample Chamber. Replaced discolored tubing and inspected 5-Way Valve. Verified Modem with laptop and DMT Tech Software using phone cord.					
Parts:	Otv	Part Number	Part Description			
	1	90-8000-00	CONTRACT PREV MAINT DMT			
Reason:	n: DMT Status Messages					
	CommEr					
Reason Notes:	son Notes: Communications Error					
Resolution:	Replaced					
	RPL comp					
Resolution Notes:	Replaced bad U31 on SI Board.					
	2	59-0085-00	DMT HRDW BASE RUBBER BUMPER 31360			
Reason:	n: DMT Functional					
	DigPots					
Reason Notes:	Requires TB-30 Digital Pots Mod.					
Resolution:	Update					
	DMTVR3					
Resolution Notes:	Removed Control Board and installed TB-30 DMT Digital Pots Mod. Optimized detector position. Adjusted Tech Screen setting					
	1	59-0036-00	DMT HRDW CLIP BREATH TUBE HOLD 37058			
	1	59-0058-00	DMT ELEC IC RS232 MULTCH U31/33			
		001199-00	DMT ELEC REGULATOR V LDO 3.3V OUT 30V IN			

## Order Packing Slip Report

Intoximeters

Experience · Service · Integrity

Pack Date: 5/26/2021

5/26/2021 10:26:33 AM

Packing Slip: 228953

SRO: SRO-046012

From: Intoximeters, Inc

2081 Craig Road St Louis MO 63146

USA

Customer: C000003651

Ship To: 2

Alaska Department of Public Safety

4805 Dr MLK Jr Ave Anchorage AK 99507

**USA** 

Alaska Department of Public Safety Procurement & Supply Section 4805 Dr. MLK Jr. Avenue Anchorage AK 99507

USA

Customer PO: 050521-Walton

Ship Via: FDX

FED EX

Weight: 37.00

Delivery Terms: 2D

FOB: ORIGIN

# Packages: 1

Customer Ship Account:

Line Number Item		UM	Ordered Quantity	Shipment Quantity	Backorder Quantity	
1	LEGDMT	LEGACY	DMT		50100352	
1	90-8000-00	EA		1.00	0.00	
CONTRACT I	PREV MAINT	DMT				
2	59-0085-00	EA		2.00	0.00	1705
DMT HRDW	BASE RUBBE	R BUMPER	31360			frig
3	59-0036-00	EA		1.00	0.00	Pec'a 05/28/2
DMT HRDW	CLIP BREATI	H TUBE HO	LD 37058			05/2010
4	59-0058-00	EA		1.00	0.00	JW
DMT ELEC IO	C RS232 MUL	TCH U31/3	3			0"
5	001199-00	EA		1.00	0.00	
DMT ELEC R	EGULATOR	V LDO 3.3\	OUT 30V IN			