Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

	ssued: 1/13/2017 Effective: 1/17/2017	Version: DMCIF 2017 R0 Status: Active		
Instru	rument S/N/00402	Received From Intoximeters		
Date_	09/15/21	Analyst Deren Watton		
1.	Reason for instrument return: Refurbished by Man	sufacturer, regard for fatal errors		
2.	. Change instrument and case location in Jus	ticeTrax Breath Alcohol Database		
3.	Ensure change in instrument status form is in	n JusticeTrax		
4.	. Upload 🗆 NA In	strument was not in the field DW 09/15/2/		
5.	. Update agency information and polling statu	s in DM Host Calling Program		
6.	. Check for VOCs			
	From Intoximeters			
7.	Intoximeters Repair Paperwork (attached)			

ntoximeters

Experience · Service · Integrity

Service Record

SRO Number: SRO-047140

Received Date: 07/26/2021

Completion Date: 08/12/2021

SRO Type: REPAIR

SRO Description: DMT Repair/Preventative Maintenance

Customer ID: C000003651

Customer PO: 07262021-Lee

Contact: Tatiana Lee

Email: tatiana.lee@alaska.gov

Ship To:

Units on SRO

Alaska Department of Public Safety Procurement & Supply Section

LEGDMT

LEGACY DMT

4805 Dr. MLK Jr. Avenue Anchorage AK 99507

Service Line:	1		**All Instruments are Calibrated to Factory Specifications**			
Unit Type:	(LEGDMT) LEGACY DMT					
Serial Number:	50100402					
Reason	DMT Functional					
	Lock Up					
Reason Notes:	Customer	Reported Issue: Locks Up or Fr	reezes, Wont Save Test Results			
Resolution	n: Adjusted					
	Mech adj					
Resolution Notes:	Added 1/16" spacer between touch screen and front panel					
Parts	Qty	Part Number	Part Description			
	1	59-0097-00	DMT HRDW SHIM WASHER SMALL 30487			
Reason: Preventive Maintenance						
Reason Notes:	rason Notes: 90-8000-00 Prev Maint Contract					
		ev M				
Resolution:	Prev M					

Resolution Notes:

90-8000-00 Prev Maint Contract p/n

Installed new BT o ring

Tightened rear panel mount connectors

Installed TB-30 modifications to the control board. SI green board B1 battery voltage checked = 3.28vdc Replaced standard slotted screws on the SC with torx screws Installed metal seal tape around the SC mirror blocks

Cleaned and adjusted the 5-way valve.

Cleaned, lubricated, & adjusted filter wheel moving parts

Adjusted the air volume: 3.01 Adjusted tech screen settings Adjusted RFI sensitivity Tested modem line pick up Time and date verified. Dry Gas leak checked

(1) 59-0097-00 Shim Washer small (1) 001199-00 3.3V regulator

Calibration solution (0.100) Lot# 21080 Bottle# 618 Dry gas: AG106803, Tank 001, Exp. 09-Mar-2023

1	001199-00	DMT ELEC REGULATOR V LDO 3.3V OUT 30V IN		
1	90-8000-00	CONTRACT PREV MAINT DMT		

Order Packing Slip Report

Intoximeters

Experience · Service · Integrity

Pack Date: 8/18/2021

8/18/2021 9:54:38 AM

Packing Slip: 234771

SRO: SRO-047140

From: Intoximeters, Inc

2081 Craig Road St Louis MO 63146

USA

Customer: C000003651

Ship To: 2

Alaska Department of Public Safety

4805 Dr MLK Jr Ave Anchorage AK 99507

USA

Alaska Department of Public Safety Procurement & Supply Section 4805 Dr. MLK Jr. Avenue Anchorage AK 99507

USA

Customer PO: 07262021-Lee

Ship Via: FDX

FED EX

Weight: 38.00

Delivery Terms: GND

Shipping Terms: ORIGIN

Packages: 1

Customer Ship Account:

Line Number	Item	UM	Ordered Quantity	Shipment Quantity	Backorder Quantity	
1	LEGDMT	LEGACY	DMT		50100402	
1	59-0097-00	EA		1.00	0.00	
DMT HRDW	SHIM WASHE	R SMALL	30487			
2	001199-00	EA		1.00	0.00	
DMT ELEC R	EGULATOR V	/ LDO 3.3V	OUT 30V IN			
3	90-8000-00	EA		1.00	0.00	
CONTRACT I	PREV MAINT	DMT				Reid 08/20/21

JM