## Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form

Issued: 7/10/2015

Effective: 7/10/2015

Instrument S/N / 100687

Date 10/15/15

Received From 1/PAS

Analyst College OBryant

1. Reason for instrument return:

Undefined error and other software (sover)

- 2. Change instrument and case location in JusticeTrax Breath Alcohol Database
- 3. Ensure change in instrument status form is in JusticeTrax
- 4. Upload MIA
- 5. Update agency information and polling status in DM Host Calling Program A/A
- 6. Check for VOCs MIA

## From NPAS 🗹

7. NPAS Repair Paperwork (attached)

Vork Ord No 49925  Date 9/4/15		National Patent Analytical Systems, Inc. Service Action Report		Serial N	Serial No 100687	
Customer	Alaska Dept of Pub	olic Phone 907	'-269 <b>-</b> 5582 <sup>Co</sup>	ntact Bolz		
War Exp	10/30/13	Мх Ехр	Ship Damage?	Y 🛈 N Cust Pick U	<sup>977</sup> : ○ Y <b>⊙</b> N	
shipping dan	ents: note any nage or estimates ng with accessories					
Describe problem seen and/or attach documents received		fatal errors afte	r software update			
Keyword		•	•			
Finding	*	Error Code	Enter Hour	s 2.0		
Repair Comr Performed	nents and work U	pdated to current terference on dry	software. Replaced de gas. Calibrated, ran s	efective touchscreer tandard tests ok	, blank errors/	
Solution Lo	t No 15040	i				
1	Voltage Set 1	Values in	Voltage Set Valu	es Out		
	Lamp 1.640	8rth Vol 1.470	Lamp 1.670	Brth Vol 1.480		
	Meter 1.530	Filter Dif 0.768	Meter -0.057			
3	Stability 0.002		Stability 0.002	Filter Dif 0.812		
Part No	Qty Ord Pick	Back Total Ordered shin	Description		Price Ea	
37553	1,0	WARR	MEMBRANE,8.4 IN TOUCH ANTY	REPLACEMENT	0.000	
9DMT	1.0	10068	DATAMASTER DMT SERVIORS	CE RETURN	0.000	
NOTE	1.0		NANUK CASE		<del> </del>	
R500-1	1.0	And the state of t	Service Charge for Repair		0.000	
R500-2	1.0		Service Charge for Repair,	Additional Hours	0.000	
<u> </u>						

Pick up Signature