

# Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form

Effective: 3/20/2025

Version: 4.0

Instrument #: 100390

Received From: Suneau PD

Date: 3/23/26

Analyst: Kaya Carlson

Reason for instrument return:

Filter wheel errors.

### From the field

Change instrument and case location in JusticeTrax Breath Alcohol Database

Ensure change in instrument status form is in JusticeTrax

Upload

Update agency information and polling status in DMHost Calling Program

Check for VOCs

### From Intoximeters

Intoximeters Repair Paperwork (attached)

Change instrument and case location in JusticeTrax Breath Alcohol Database

Reset Options and Clear Memory