VERIFICATION OF CALIBRATION REPORT

of DataMaster dmt Breath Test Instrument State of Alaska

Scientific Crime Detection Laboratory - Statewide Breath Alcohol Program

Date:01/24/2013

External Standard Test Values

EXTERNAL STANDARD INFORMATION

NOMINAL: 0.080

Serial #: 100420

TARGET AT 28.32: 0.075 LOT #: 22612080A1 EXPIRATION: 09/01/2014

BLANK TEST	0.000	18:28
INTERNAL STANDARD	VERIFIED	18:28
EXTERNAL STANDARD	0.075	18:28
BLANK TEST	0.000	18:30
EXTERNAL STANDARD	0.076	18:30
BLANK TEST	0.000	18:31
EXTERNAL STANDARD	0.076	18:31
BLANK TEST	0.000	18:32
EXTERNAL STANDARD	0.076	18:32
BLANK TEST	0.000	18:33
EXTERNAL STANDARD	0.076	18:33
BLANK TEST	0.000	18:34

Average = 0.0758Std Dev = 0.0004

Diagnostic Check

VERSIONS DMT: 2.00 PIC: 3.02 Modem: 2.3 Questions: 2.2

INTERNAL STANDARD

TEMPERATURES	
Sample Chamber = 49.0°C	PASSED
Breath Tube = 43.4°C	PASSED
PUMP INFO	
Flow Rate = 4.965 L/M	PASSED
DETECTOR INFO	
PUMP ON	PASSED
PUMP OFF	PASSED
FILTER INFO	
Filter 1	PASSED
Filter 2	PASSED
Filter 3	PASSED

I, Nita J. Bolz, after being first duly sworn, depose and state as follows:

(1) I am a Forensic Scientist IV at the State of Alaska Scientific Crime Detection Laboratory.

(2) The Alaska Scientific Crime Detection Laboratory is an entity within the Department of Public Safety.

(3) I am the Scientific Director of the State Breath Alcohol Program.

(4) In that capacity, I am responsible for overseeing the Breath Alcohol Program, which includes assuring that instruments are calibrated and maintaining program records.

(5) The above is a true and accurate verification of calibration, which is performed by the instrument's software, as specified by the State Breath Alcohol Program. Verification of calibration is a regularly conducted and regularly recorded activity of the State Breath Alcohol Program.

(6) The referenced instrument is certified for evidentiary use in the State of Alaska.

Nita J.Bolz

Scientific Director State Breath Alcohol Program

Subscribed and sworn before me this

Nikki Roth, Notary Public My Commission Expires With Office





PASSED