

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Effective: 10/10/2022

Version: 3.0

Instrument #: 100W07

Received From: North Pole PD

Date: ~~01/06/22~~ 12/06/22
FLM 12/07/22

Analyst: Tatiana Maloske
Supervised by: Colleen O'Bryant

Reason for instrument return:

Accidentally pulled for PM. Instrument only in service
for 5 months. calibrated instrument upon removal.

Change instrument and case location in JusticeTrax Breath Alcohol Database

Ensure change in instrument status form is in JusticeTrax

Upload

Update agency information and polling status in DMHost Calling Program

Check for VOCs

From Intoximeters

Intoximeters Repair Paperwork (attached)

*Tightened top right
screw for touch screen
to stop triggering the
tank change button
w/out pressing.
C.03 12/6/22*