## Alaska Scientific Crime Detection Laboratory

## DataMaster Evaluation Form Issued: 05/15/2020 Version: DMEF 2020 R0 Effective: 05/15/2020 Status: Active Instrument # /00352 Analysts: Derch walter Reason for instrument evaluation Preventive maintenance duc Visual Inspection Replace missing screws Clean Cermetex Box is secure Ensure modem cable present Gas Valve upright or secured to wall All tubing and connections are secure O-ring inspection and ensure 2 O-rings in replacement bag Troubleshooting/Repair Added moder cable, Replaced o-ring. Tightened source block. Now snowing communication ecross. Instrument also see needs VOC (attached) No further tisting Im 05/23/21 Internal Standard value ≤ 3.00% Yerify software version Return to Factory for Repair Ready for Calibration Ready for Certification Proceed to Close Out End Date: 03/23/21 Page 1 of 1

Approved by Chemistry Supervisor

All printed copies are uncontrolled.