

Alaska Scientific Crime Detection Laboratory

DataMaster Check In Form

Effective: 3/20/2025

Version: 4.0

Instrument #: 100365

Received From: Eielson AFB

Date: 06/26/2025

Analyst: Dereen Watten

Reason for instrument return:

Instrument never placed in service after calibration due to
Interference Detected and Internal Standard errors on
VOCs

From the field

Change instrument and case location in JusticeTrax Breath Alcohol Database ☒

Ensure change in instrument status form is in JusticeTrax ☒

Upload ☐ Could not upload via modem line. Output files to USB JMW 06/26/25

Update agency information and polling status in DMHost Calling Program ☒

Check for VOCs ☒

From Intoximeters

Intoximeters Repair Paperwork (attached) ☐

Change instrument and case location in JusticeTrax Breath Alcohol Database ☐

Reset Options and Clear Memory ☐