## Alaska Scientific Crime Detection Laboratory

## DataMaster Check In Form

Issued: 1/13/2017 Effective: 1/17/2017

Version: DMCIF 2017 R0

Status: Active

Instr	ument S/N 100361 Received From Unalaska AD
Date_	
1.	Reason for instrument return:  Preventative Maintenance
2.	Change instrument and case location in JusticeTrax Breath Alcohol Database
3.	Ensure change in instrument status form is in JusticeTrax
4.	Upload
5.	Update agency information and polling status in DM Host Calling Program
6.	Check for VOCs

## From Intoximeters

7. Intoximeters Repair Paperwork (attached)