Alaska Scientific Crime Detection Laboratory

DataMas	ster Evaluation Form
Issued: 05/15/2020 Effective: 05/15/2020	Version: DMEF 2020 R0 Status: Active
Instrument #/00699	
Start Date:03/08/21	Analysts: Derch Walton
Reason for instrument evaluation Blank errors on	attempted calibration (H20)
Visual Inspection Clean	□ Replace missing screws Sel evaluation from
☐ Cermetex Box is secure	Ensure modem cable present /2/18/20 m/
☐ Gas Valve upright or secured to wa	II \square All tubing and connections are secure 03/18/21
O-ring inspection and ensure 2 O-ri	ngs in replacement bag
No Further translesshood	ting performed factory
VOC (attached) N → perfilms Internal Standard value ≤ 3.00% Verify software version	ed 03/08/21
Return to Factory for Repair	☐ Ready for Calibration
Ready for Certification	Proceed to Close Out
End Date: 03/08/21	
All printed copies are uncontrolled.	Page 1 of 1 Approved by Chemistry Supervisor