Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form Issued: 04/03/2018 Version: DMEF2018 R0 Effective: 04/04/2018 Status: Active Instrument #_____100342 Date: 09/08/20 Analyst: Durch Watton 1. Reason for instrument evaluation Static analysis errors 2. Visual Inspection Clean Replace missing screws ☐ Cermetex Box is secure Ensure modem cable present ☐ Gas Valve upright or secured to wall ☐ All tubing and connections are secure 3. O-ring inspection and ensure 2 O-rings in replacement bag 4. Troubleshooting/Repair 5. VOC (attached) 6. Internal Standard value ≤ 3.00% Yes □ No \square 7. Verify software version Return to Factory for Repair Ready for Calibration Ready for Certification ☐ Proceed to Close Out