Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form Effective: 10/10/2022 Version: 3.0 Instrument #: __/0069 Start Date: __ \Q/14/6 Date Received: Reason for instrument evaluation Touch screen not responding Visual Inspection Replace missing screws Clean Cermetex Box (modem) is secure ☑ Ensure modem cable present ☐ Double Check Valve secured to wall ☑ All tubing and connections secure O-ring inspection and ensure 2 replacement O-rings in bag As Found Conditions: Temperature (C):_____Relative Humidity (%)_____ Barometric Pressure (inHg): COB 10/14/00 Linearity Test (attached) 0.020g/210L Mnf: Lot: Lot: 0.080g/210L Mnf:___ 0.300g/210L Mnf:_ Lot: Notes unable to perform au-found due to lack of VOC (attached) Verify software version VOC Internal Standard value ≤ 3.00% Non-Drinking Subject Test (attached) Needs Repair Condition of item: Good Ready for Adjustment ☐ Ready for Calibration

Page 1 of 1

All printed copies are uncontrolled.

Approved by: Chemistry Supervisor