Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form Issued: 04/03/2018 Version: DMEF2018 R0 Effective: 04/04/2018 Status: Active Instrument # 100c I. Reason for instrument evaluation Pump failure on diagnostics + fatal system 2. Visual Inspection ☐ Clean Replace missing screws ☐ Cermetex Box is secure Ensure modem cable present ☐ Gas Valve upright or secured to wall All tubing and connections are secure 3. O-ring inspection and ensure 2 O-rings in replacement bag 4 Troubleshooting/Repair The screen is having issues and only shows a Recommending refurbishment @ Factory 5. VOC (attached) 6. Internal Standard value ≤ 3.00% Yes □ No 🗔 7. Verify software version Return to Factory for Repair Ready for Calibration Proceed to Close Out Ready for Certification