Alaska Scientific Crime Detection Laboratory

DataMaster Evaluation Form Issued: 04/03/2018 Version: DMEF2018 R0 Effective: 04/04/2018 Status: Active Instrument # 106691 Analyst: Dereh Walton 1. Reason for instrument evaluation Frequent freezing, updated software did not help 2. Visual Inspection ☐ Clean Replace missing screws Ensure modern cable present ☐ Cermetex Box is secure All tubing and connections are secure ☐ Gas Valve upright or secured to wall 3. O-ring inspection and ensure 2 O-rings in replacement bag 4. Troubleshooting/Repair 5. VOØ (attached) 6. Internal Standard value ≤ 3.00% Yes No \square 7. Verify software version Return to Factory for Repair Ready for Calibration Proceed to Close Out Ready for Certification