## Alaska Scientific Crime Detection Laboratory Change in Instrument Status Form

Effective: 8/10/2022

Version: 4.0

tersioni no	2110001101 0/ 20/ 2022
Date: 2/17/2023	Instrument #: 100403
nstrument Location: Wasilla PD	
Supervisor Name and Agency: Colleen O'Bryant, SCDL	
Describe the Status Change or Issue with t	he Instrument:
· · · · · · · · · · · · · · · · · · ·	or status message at Wasilla PD instrument. He was able to data screen in a subject test, selecting the drop down for test type pace.
This is a software issue the Breath Program i type from the drop down rather than draggi	is aware of. The status message can be avoided by selecting a test ng the stylus down the options.
If instrument was placed in service or removed from service, fill out the appropriate sections below.	
In Service Date/Time: N/A	
Out of Service Date/Time: N/A	
For Use by SCDL	Additional Notes
Click or tap here to enter text.	

**Email completed form to** <u>dps.scdl.toxicology@alaska.gov</u>
For questions contact the Breath Alcohol Program at 907-269-5740

Page 1 of 1

All printed copies are uncontrolled.

Approved by: Chemistry Supervisor

## O'Bryant, Colleen S (DPS)

From: Adams, Andrew J (DPS)

Sent: Friday, February 17, 2023 8:45 PM

**To:** O'Bryant, Colleen S (DPS) **Subject:** Re: Datamaster error WPD

Very interesting.

Got it to do it on another Datamaster.

Neat party trick.

On Feb 17, 2023, at 19:03, O'Bryant, Colleen S (DPS) <colleen.obryant@alaska.gov> wrote:

Start a test. From the data entry screen, select the drop down for test type and hover your finger/stylus on the blank space.

## Colleen O'Bryant

Forensic Scientist III - Chemistry Alaska Scientific Crime Detection Laboratory 4805 Dr. Martin Luther King Jr. Ave Anchorage, Alaska 99507

907-269-5783 (direct) 907-350-2111 (cell)

From: Adams, Andrew J (DPS) <andrew.adams@alaska.gov>

**Sent:** Friday, February 17, 2023 7:02:15 PM

To: O'Bryant, Colleen S (DPS) <colleen.obryant@alaska.gov>

Subject: Re: Datamaster error WPD

Interesting. I might try and recreate that on a different Datamaster.

Sent from my iPhone

On Feb 17, 2023, at 17:21, O'Bryant, Colleen S (DPS) <colleen.obryant@alaska.gov> wrote:

I have. It usually happens when you hover too long on the test type.

## Colleen O'Bryant

Forensic Scientist III - Chemistry Alaska Scientific Crime Detection Laboratory 4805 Dr. Martin Luther King Jr. Ave Anchorage, Alaska 99507

```
907-269-5783 (direct)
907-350-2111 (cell)
```

From: Adams, Andrew J (DPS) <andrew.adams@alaska.gov>

Sent: Friday, February 17, 2023 4:52:12 PM

To: O'Bryant, Colleen S (DPS) <colleen.obryant@alaska.gov>

Subject: Datamaster error WPD

I have never seen this error before.

<image0.jpeg> <image1.jpeg>